



## DynaPlus™ Nitrile Exam Gloves




## dynap<sup>+</sup> Nitrile Gloves Provide Strength Plus Value!

*DynaPlus™ Nitrile Exam Gloves deliver the protection you need in a value pack –Double the quantity in box. Value pack saves you time and money!*

- Superior Comfort • Textured fingertips
- Strong, easy-to-don beaded cuffs
- Powder-Free • Latex Free

Reorder No.	Description	Quantity
2516	Nitrile Exam Glove - S	10/200/Cs
2517	Nitrile Exam Glove - M	10/200/Cs
2518	Nitrile Exam Glove - L	10/200/Cs
2519	Nitrile Exam Glove - XL	10/180/Cs

Color	Thickness	Textured	Tensile Strength	Cuff	Meets Standards	Market Position
	3 mil	Fingertip	14 MPa, Min.	Beaded	ASTM D6319, ASTM D6124	Economy



Already mobile?  
Scan your smart phone here for more information.

See our website for our High Risk Latex Gloves (2310-2318)

**Find us on:**



For more information contact your local Dynarex Distributer Sales Representative.  
Or call 888-396-2739 to find your nearest distributor.