



WE CARE LIKE FAMILY



Compares to the Welch Allyn Universal Kleenspec Disposable Otoscope Specula #52432 & 52434*

Disposable Ear Specula Tips

Fits all major brand diagnostic otoscopes on the market

- *High quality medical grade plastic with a secure fit design and no sharp edges*
- *Disposable tips*
- *Economical and convenient*
- *Not made with natural rubber latex*

7116
2.5 mm



7117
4.25 mm



NOTE: Dynarex Otoscope tips fit the following Welch Allyn Diagnostic Otoscopes: Item #20000, 25020, 24031. They do not fit MacroView™ Oscopes.

Reorder No.	Description	Quantity
7116	2.5 mm, Pediatric	10/850/Cs
7117	4.25 mm, Adult	10/850/Cs

*Welch Allyn Universal Kleenspec Disposable Otoscope Specula and MacroView are trademarks of Welch Allyn.



See our website for our Portable Penlights (#4781-4782)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com